Laborporzellan / *Laboratory Porcelain*

ALSINT 99,7 Oxidkeramik / Oxide Ceramics

CaO-FSZ Zirkonoxid / Zirconia

Inhalt / Content

Laborporzellan	Seite	Laboratory Porcelain	Page
Werkstoffe	5	Materials	5
Mörser 55	10	Mortars 55	10
Pistillen 56	10	Pestles 56	10
Abdampfschalen 8, 8 A	10	Evaporating basins 8, 8 A	10
Abdampfschalen 109, 130, 131, 888	11	Evaporating basins 109, 130, 131, 888	11
Glühschalen 33, 33 C, 33 D	12	Incinerating dishes 33, 33 C, 33 D	12
Mehlveraschungsschalen 5032	12	Flour incinerating dishes 5032	12
Schmelztiegel 79, 79 C, 79 MF	13	Melting crucibles 79, 79 C, 79 MF	13
Filtriertiegel 84	13	Filtering crucibles 84	13
Gooch-Tiegel 82 A, 82 R	14	Gooch crucibles 82 A, 82 R	14
Siebplatten 31 B	14	Filtering discs 31 B	14
Tiegeldeckel 79 D	14	Crucible lids 79 D	14
Platten 61 A, 61 B	15	Spot plates 61 A, 61 B	15
Exsikkatoren-Platten 119 C	15	Desiccator plates 119 C	15
Kasserollen 17	15	Casseroles 17	15
Trichter 126, 127 C	16	Funnels 126, 127 C	16
Schiffchen 30	16	Combustion boats 30	16
Messbecher 51	17	Graduated beakers 51	17
Kochbecher 42 A, 42 L	17	Heating beakers 42 A, 42 L	17
Kugelmühlen GSK und Kugeln 42 K	18	Ball mills GSK and balls 42 K	18
Löffelspatel 74	19	Spoon spatulas 74	19
Tonteller	19	Clay plates	19
Hartporzellan-Rohre	19	Hard porcelain tubes	19
Porzellan-Schiffchen	19	Porcelain combustion boats	19
Alsint 99,7 Oxidkeramik		Alsint 99.7 Oxide Ceramics	
Werkstoffe	5	Materials	5
Alsint 99,7-Schiffchen	20	Alsint 99.7 Boats	20
Alsint 99,7 porös-Glühschalen	20	Alsint 99.7 porous Incinerating dishes	20
Alsint 99,7-Glühkästen	20	Alsint 99.7 Incinerating dishes	20
Alsint 99,7-Röhrentiegel	20	Alsint 99.7 Tubular crucibles	20
Alsint 99,7-Tiegel	21	Alsint 99.7 Crucibles	21
Zirkonoxid CaO-FSZ		Zirconia CaO-FSZ	
Werkstoffe	5	Materials	5
CaO-FSZ-Tiegel	22	CaO-FSZ Crucibles	22

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Haldenwanger Laboratory Porcelain, the laboratory specialist – since 1865

We have risen to many different challenges in over 150 years:

High Temperatures

	Heat resistant up to max. temp.:
Glazed laboratory porcelain:	1,000 °C
Unglazed laboratory porcelain:	1,350 °C
Alsint 99.7:	1,800 °C
CaO-FSZ:	2,000 °C

Aggressive Reagents

Laboratory Porcelain	Chemical resistance:
All acids (20 °C):	excellent*
All acids (boiling):	excellent*
Warm, highly-concentrated alkaline solutions:	good
Molten alkaline salts:	unsuitable
	*Except hydrofluoric acid

Alsint 99.7

Reliable resistance to diverse, in some cases highly-concentrated reagents both at room temperature and boiling point, except the following: highly concentrated hydrofluoric acid; boiling phosphoric acid; boiling potassium hydroxide and sodium hydroxide solutions.

Zirconia CaO-FSZ

High refractory properties in temperatures of up to 2,000 °C. In comparison with aluminium oxide improved chemical resistance to alkalis, acids and caustics.

High-wear Conditions

For example: mortar and pestle, ball mill, mouthpieces, shaft-protection sleeves, etc.

New Developments

We are constantly optimising our laboratory porcelain to meet your needs.

Highest Standards

Our laboratory porcelain complies with DIN EN 60672-3, Group C100, Type C110; laboratory articles made of Alsint 99.7 comply with DIN EN 60672, Group C700, Type C799. Stringent product controls, from raw materials to the finished product, ensure the consistently high quality of our products.

You will find applications and directions for use on our website: <http://www.haldenwanger.de/>

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Please note that all values quoted are based on test specimens and may vary according to component design. These values cannot be guaranteed and can only be transferred to other forms and dimensions to a limited extent. They should be used for guidance only. In the field, for example, Alsint 99.7 moulded pieces demonstrate mechanical strength values between 160 and 300 MPa, depending on wall thickness, actual shape, surface finish, the shaping process and post-processing.

* We recommend that product be heated at a rate not exceeding 30-50 °C per hour
** dependent on load

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
Why should I use Haldenwanger laboratory porcelain?

Haldenwanger has been a producer of laboratory porcelain for 150 years.

Despite such a long history record, we have taken care to maintain the original size and shape of our products. New challenges of the ever growing market have led us to new developments. The best refractoriness and resistance against physical and chemical attack has always been of utmost relevance to us.

Technical Support

As a producer of high performance ceramics, Haldenwanger offers not only the best materials but also technical support on materials selection and application. In the area of laboratory ceramics, we have a product portfolio ranging up to 2,000°C temperature stability.

Werkstoff / Material	Anwendungstemperatur / Application temperature	WAK*10 ⁻⁶ K ⁻¹
Porzellan glasiert / Porcelain glazed	1000°C	5,3
Porzellan unglasiert / Porcelain unglazed	1350°C	5,3
Al ₂ O ₃	1800°C	8-9
ZrO ₂ -CaO-FSZ	2000°C	10

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Being a member of the VGKL (Verband des Groß- und Außenhandels für Krankenpflege- und Laborbedarf), a trade association of the leading wholesalers for laboratory equipment, we have always been at the cutting edge of technology. In addition, Haldenwanger offers training courses in the theory and usage of laboratory porcelain.

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Proven Quality of Haldenwanger Ceramics

We guarantee a consistently high quality level by performing quality checks at all levels, e.g. raw materials, forming process, and final product inspection. To this end, we have established our own state-of-the-art testing and development centre as well as a pilot plant. Our products are continuously tested to comply with DIN 12851 regulations.

This DIN standard comprises substantial and strict requirements for the testing of laboratory equipment made from alumo-silicate ceramics (materials class C110). In particular, this encompasses the following requirements:

- *imperviousness of the ceramic body*
- *dense glazing (no pores, no cracks)*
- *thermal shock resistance*
- *refractoriness of the glazing up to 900 °C*
- *no weight change during annealing*
- *chemical resistance against hydrochloric acid, diluted soda and caustic soda*

Requirements

Imperviousness of the ceramic body

No single area of any sample may show

- *any kind of discoloration in glaze-free areas*
- *penetration of dye between ceramic body and glazing*

Dense glazing (no pores, no cracks)

No single sample may show trajectories of excess glazing slurry or any discoloration along the edges.

Thermal shock resistance

No single sample may fail by fracture and no crazing may appear after quenching.

Refractoriness of the glazing up to 900 °C

Glazed pieces may not stick together at 900 °C.

No weight change during annealing

No single sample may show a weight change exceeding 0.1 mg per 10 g material.

Chemical resistance against hydrochloric acid, soda and caustic soda

No single sample may show a weight loss

- *of more than 1 mg per 1 dm² of the wetted inner surface area (hydrochloric acid)*
- *of more than 10 mg per 1 dm² of the wetted inner surface area (diluted soda)*
- *of more than 60 mg per 1 dm² of the wetted inner surface area (caustic soda)*

Haldenwanger Quality Proven to be the Best

The following images show a Haldenwanger as well as a competitor's product before and after single usage at 1,000 °C.

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Haldenwanger crucibles shown before (left) and after (right) anneal testing. The result confirms the thermal stability of the glazing, even after multiple testing cycles.

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Competitor crucibles shown before (left) and after (right) anneal testing. Significant discolorations can be detected after the first anneal cycle.

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
Why is thermal shock resistance of utmost relevance for laboratory porcelain?

Thermal shock resistance as described in the DIN 51068 is the material's response to repeated subjection to thermal stresses occurring in (rapidly) changing temperature fields. Heating of laboratory ceramics and cooling may lead to thermal stresses within the ceramic body. Due to the temperature gradients within the material, compressive and tensile stresses may form. In extreme situations, local stresses exceed the strength level of the material, resulting in cracking.

Thus excellent thermal shock properties are most relevant for longevity of the laboratory porcelain.

The following table shows the result of a thermal shock test performed at 1,000 °C on a Haldenwanger and several competitors' products.

Test product: is equivalent to Haldenwanger crucible type 79 MF/6

Test conditions: empty crucibles were placed in a pre-heated furnace at 1,000 °C, anneal time was 15 min, then removal from the furnace and cool-down

Testing: running five cycles of heating-up and cooling

Analysis: proportion of crucibles without fracture

In order to extend the life-time expectancy of your laboratory porcelain, we recommend you follow these instructions:

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Morgan Advanced Materials Haldenwanger is your number one specialist in high-temperature and laboratory ceramics. Having gained experience over more than 150 years, we can offer the ideal solution. We would be pleased to advise you – in many European languages!

Hints for Using Laboratory Porcelain Correctly

Life-time expectancy of laboratory porcelain is a function of various factors such as heat-up and cooling rates, degree of filling (half/completely filled), homogeneity of the filling within the crucible, geometry of the crucible (round / sharp curved, wall thickness, radii, wall thickness transitions), usage (e.g. during cleaning), etc. In order to avoid damage of the material, the user should respect the following guidelines:

Maximum temperature

- Glazed porcelain: 1,000 °C
- Unglazed porcelain: 1,350 °C
- Al₂O₃ crucibles: 1,800 °C
- CaO-FSZ crucibles: 2,000 °C
- Maximum heating rate of porcelain crucibles: 150 °C/h
- Maximum heating rate of Al₂O₃ and CaO-FSZ crucibles: 50 °C/h
- At higher application temperatures (> 400 °C), the curved crucibles 33 and 33D should preferably be used over the sharp-edged type 33C
- Choosing the right size compared with the content, minimum degree of filling should exceed 75% (to avoid temperature differences)
- Homogeneous filling within the crucible

Operation during cool-down phase

- only use pre-heated crucible tongs
- place hot crucibles only on porous ceramic materials such as R-SiC or alumo-silicate materials

Generally, fast temperature cycling should be avoided. Most critical is heat loss during exceedingly fast cooling.

The instructions given above are based on our long-term experience with our products. Further factors such as furnace atmosphere or chemical interactions with the material treated may also influence the thermal shock resistance of porcelain.

Custom-made Ceramics

The following pages will give you an overview of the current product portfolio of our laboratory porcelain. In addition, we offer individualised products in variable material and shape/size combinations. Also, we can apply your logo, any serial number, or any bar codes on the glazing.

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Cylindrical, with shallow base, with or without spout, glazed except outside base
Without spout also usable as lid for article L42 A or L42 L

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Form B, semi-deep, with spout, Size No. 5/0 – 5 glazed inside and outside, Size No. 6 – 10 glazed except outside base

Größen-Nr. Size No.	DIN	Ø mm	Höhe / Height mm	Inhalt / Capacity ca. / approx. ml	VE Sales Unit
5/0	I 2903	40	16	10	100
4/0	I 2903	50	20	20	100
000	I 2903	63	25	40	100
00		70	26	50	100
0	I 2903	80	34	75	100
1		85	35	95	100
2	I 2903	100	40	115	50
3		105	42	170	50
4		110	45	220	50
5	I 2903	115	47	310	10
6		150	45	370	10
6 a	I 2903	160	64	640	5
7 b	I 2903	200	80	1000	5
9	I 2903	200	86	2200	2
10	I 2903	310	110	3400	2

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french shape, with spout and round base, glazed inside and outside

Größen-Nr. Size No.	Ø mm	Höhe / Height mm	Inhalt / Capacity ca. / approx. ml	VE Sales Unit
3	75	40	70	10
5	97	50	190	10
6	110	52	270	10
7	125	60	350	10

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french shape, with spout and flat base, glazed inside and outside

Größen-Nr. Size No.	Ø mm	Höhe / Height mm	Inhalt / Capacity ca. / approx. ml	VE Sales Unit
2	58	24	30	20
3	75	30	60	20
4	84	36	100	10
5	97	40	150	10
6	110	50	250	10
7	125	52	290	10
8	150	60	400	5

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Form A, shallow, with spout, glazed except outside base

Größen-Nr. Size No.	DIN	Ø mm	Höhe / Height mm	Inhalt / Capacity ca. / approx. ml	VE Sales Unit
000	I 2903	40	9	5	30
00	I 2903	50	11	10	30
0	I 2903	63	13	20	30
2	I 2903	80	20	40	20
4	I 2903	100	22	80	20
6	I 2903	125	27	270	20
6 a	I 2903	160	35	450	10
8	I 2903	190	55	1100	3
9	I 2903	240	65	1750	3
10	I 2903	300	60	2500	2

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Ø mm	Höhe / Height mm	Inhalt ca. ml Capacity approx. ml	VE Sales Unit
56	23	30	100

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Lid 79 D, filtering discs 3l B have to be ordered separately.
Filter holes approx. 0.5 mm in Ø

Artikel Nr. / Article No. 82 R
hohe Form, mit Siebboden, glasiert
tall shape, with perforated base, glazed

Größen-Nr. Size No.	oberer / upper Ø mm	Boden / Base mm	Höhe / Height mm	Inhalt / Capacity ca. / approx. ml	VE Sales Unit
3	35	24	41	25	20

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* The nominal size is equivalent to the diameter of suitable filter paper.

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Custom Design

We can custom-make Alsint 99.7 components to your requirements. Simply send us your detailed specifications and we will provide you with an obligation-free quote. Alsint 99.7 is the material of choice for use under high-wear conditions, in chemical technical and electrical engineering applications and in high-temperature technology.

Important Information Concerning the Use of Alsint 99.7 Laboratory Equipment:

Due to its high refractoriness (melting point above 2,000 °C) and chemical resistance to a wide range of substances, crucibles and other devices made of Alsint 99.7 are used in multiple applications. Its outstanding properties, however, are only fully maintained with proper handling.

Applications

Chemical Technical Applications:

Crucibles for crystal growing, crucibles for fusion processes, crucibles for annealing

High-wear Conditions:

Ball mill pots, milling balls, mouthpieces, shaft-protection sleeves

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Special sizes and custom designs available on request.

All products listed in this catalogue are in accordance with DIN
40680 – subject to inadvertent modifications, errors and omissions.

Morgan Advanced Materials
Haldenwanger GmbH
Teplitzer Str. 27
84478 Waldkraiburg
Germany

T +49 (0)8638 60 04-0
F +49 (0)8638 60 04-317

info.haldenwanger@morganplc.com

<http://www.haldenwanger.de>

Morgan Advanced Materials Haldenwanger ist international zertifiziert nach ISO 9001:2008, ISO 14001:2004, BS OHSAS 18001:2007
Morgan Advanced Materials Haldenwanger has achieved international ISO 9001:2008, ISO 14001:2004 and BS OHSAS 18001:2007 certification.

www.morganadvancedmaterials.com
www.haldenwanger.de